

Macconkey Agar “TBL”

Brand: Hetu

Product Code: MAA020

Availability: 10

Weight: 0.00kg

Dimensions: 0.00cm x 0.00cm x 0.00cm

Price: ?468

Ex Tax: ?468

Short Description

Macconkey Agar “TBL” *For cultivation and differentiation of enteric bacteria and gram-positive microorganisms from clinical samples especially Enterococcus faecalis within 18 hours.*

Specification

Product	
Pack Size	10 plates
Pack Size	20 plates
Pack Size	50 plates